

## ZBTB8A antibody

### Product Information

Catalog No.:	FNab09599
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in transcriptional regulation.

### Immunogen information


Immunogen:	zinc finger and BTB domain containing 8A
Synonyms:	Zinc finger and BTB domain-containing protein 8A BTB/POZ and zinc-finger domain-containing factor BTB/POZ and zinc-finger domains factor on chromosome 1 (BOZ-F1) ZBTB8A BOZF1
Observed MW:	50 kDa, 38 kDa
Uniprot ID :	Q96BR9

### Application

Reactivity:	Human
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:500-1:5000; IF: 1:20-1:200; IP: 1:200-1:1000
Image:	



IP Result of anti-ZBTB8A (IP:FNab09599, 3ug;  
Detection:FNab09599 1:300) with HEK-293 cells  
lysate 920ug.



HEK-293 cells were subjected to SDS PAGE  
followed by western blot with  
FNab09599(ZBTB8A Antibody) at dilution of  
1:1000