

## ZC3H11A antibody

### Product Information

Catalog No.:	FNab09600
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in nuclear mRNA export; probably mediated by association with the TREX complex.

### Immunogen information

Immunogen:	zinc finger CCCH-type containing 11A
Synonyms:	Zinc finger CCCH domain-containing protein 11A ZC3H11A KIAA0663 ZC3HDC11A
Observed MW:	89 kDa
Uniprot ID :	O75152

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:20-1:200
Image:	



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab09600(ZC3H11A Antibody) at dilution of 1:600