

## ZC3H7B antibody

### Product Information

Catalog No.:	FNab09604
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	zinc finger CCCH-type containing 7B
Synonyms:	Zinc finger CCCH domain-containing protein 7B Rotavirus 'X'-associated non-structural protein (RoXaN) ZC3H7B KIAA1031
Observed MW:	112 kDa
Uniprot ID :	Q9UGR2

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



HUVEC cells were subjected to SDS PAGE followed by western blot with FNab09604(ZC3H7B Antibody) at dilution of 1:600