

ZCCHC7 antibody

Product Information

Catalog No.:	FNab09612
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	zinc finger, CCHC domain containing 7
Synonyms:	Zinc finger CCHC domain-containing protein 7 TRAMP-like complex RNA-binding factor ZCCHC7 ZCCHC7
Observed MW:	63 kDa
Uniprot ID :	Q8N3Z6

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09612(ZCCHC7 Antibody) at dilution of 1:300