

## ZDHHC15 antibody

### Product Information

Catalog No.:	FNab09616
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Palmitoyltransferase specific for GAP43 and DLG4/PSD95.

### Immunogen information

Immunogen:	zinc finger, DHHC-type containing 15
Synonyms:	Palmitoyltransferase ZDHHC15 Acyltransferase ZDHHC15 Zinc finger DHHC domain-containing protein 15 (DHHC-15) ZDHHC15
Observed MW:	42 kDa
Uniprot ID :	Q96MV8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with FNab09616(ZDHHC15 Antibody) at dilution of 1:1000