

ZDHHC9 antibody

Product Information

Catalog No.:	FNab09619
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The ZDHHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS.

Immunogen information


Immunogen:	zinc finger, DHHC-type containing 9
Synonyms:	Palmitoyltransferase ZDHHC9 Zinc finger DHHC domain-containing protein 9 (DHHC-9, DHHC9) Zinc finger protein 379 Zinc finger protein 380 ZDHHC9 CXorf11 ZDHHC10 ZNF379 ZNF380
Observed MW:	38-48 kDa
Uniprot ID :	Q9Y397

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab09619(ZDHHC9 Antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab09619(ZDHHC9 Antibody) at dilution of 1:300