

ZFP36 antibody

Product Information

Catalog No.:	FNab09623
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

mRNA-binding protein involved in post-transcriptional regulation of AU-rich element(ARE)-containing mRNAs. Acts by specifically binding ARE-containing mRNAs and promoting their degradation. Recruits deadenylase CNOT7(and probably the CCR4-NOT complex) via association with CNOT1. Plays a key role in the post-transcriptional regulation of tumor necrosis factor(TNF). Plays a key role in the post-transcriptional regulation of tumor necrosis factor(TNF).


Immunogen information

Immunogen:	zinc finger protein 36, C3H type, homolog(Mouse)
Synonyms:	mRNA decay activator protein ZFP36 G0/G1 switch regulatory protein 24 Growth factor-inducible nuclear protein NUP475 Tristetraprolin Zinc finger protein 36 (Zfp-36) ZFP36 G0S24 NUP475 RNF162A TIS11A TTP
Observed MW:	42 kDa
Uniprot ID :	P26651


Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100

Image:



Immunohistochemistry of paraffin-embedded human bladder tissue slide using FNab09623(ZFP36 Antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab09623(ZFP36 antibody) at dilution of 1:300