

## ZFP36L1 antibody

### Product Information

Catalog No.:	FNab09624
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probable regulatory protein involved in regulating the response to growth factors.

### Immunogen information


Immunogen:	zinc finger protein 36, C3H type-like 1
Synonyms:	mRNA decay activator protein ZFP36L1 Butyrate response factor 1 EGF-response factor 1 (ERF-1) TPA-induced sequence 11b Zinc finger protein 36, C3H1 type-like 1 (ZFP36-like 1) ZFP36L1 BERG36 BRF1 ERF1 RNFI62B TIS11B
Observed MW:	46 kDa
Uniprot ID :	Q07352

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab09624(ZFP36L1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09624(ZFP36L1 Antibody) at dilution of 1:300