

anti- ZFP90 antibody

Product Information

Catalog No.:	FNab09627
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May function as a repressor or silencer protein, and most likely exerts its repressing activity upon zinc-dependent binding to DNA. May be involved in proper spermatogenesis by repressing the expression of genes unnecessary or incompatible with the maintenance of a haploid cell state(By similarity). Isoform 2 acts as a bridge between FOXP3 and the corepressor TRIM28, and is required for the transcriptional repressor activity of FOXP3 in regulatory T-cells(Treg)(PubMed:23543754).

Immunogen information

Immunogen:	zinc finger protein 90 homolog(Mouse)
Synonyms:	KIAA1954, NK10, Zfp 90, ZFP90, Zinc finger protein 90 homolog, ZNF756
Observed MW:	72 kDa
Uniprot ID :	Q8TF47

Application

Reactivity:	Human, Rat, Mouse
-------------	-------------------

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889

www.fn-test.com

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab09627(ZFP90 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com