

## anti- ZGLP1 antibody

### Product Information

Catalog No.:	FNab09633
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Transcriptional repressor that plays a central role in somatic cells of the gonad and is required for germ cell development. Able to repress GATA transcription factor function(By similarity).

### Immunogen information

Immunogen:	GATA like protein-1
Synonyms:	None
Observed MW:	30 kDa
Uniprot ID :	P0C6A0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human ovary using FNab09633(GLP1 antibody) at dilution of 1:50



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab09633(ZGLP1 Antibody) at dilution of 1:300