

anti- ZHX2 antibody

Product Information

Catalog No.:	FNab09636
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as a transcriptional repressor. Represses the promoter activity of the CDC25C gene stimulated by NFYA.

Immunogen information

Immunogen:	zinc fingers and homeoboxes 2
Synonyms:	AFP regulator 1, AFR1, Alpha fetoprotein regulator 1, KIAA0854, RAF, Regulator of AFP, ZHX2, zinc fingers and homeoboxes 2
Observed MW:	92 kDa
Uniprot ID :	Q9Y6X8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:1000-1:10000; IP: 1:200-1:2000
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP Result of anti-ZHX2 (IP:FNab09636, 4ug;
Detection:FNab09636 1:800) with mouse brain
tissue lysate 2800ug.



human heart tissue were subjected to SDS PAGE
followed by western blot with FNab09636(ZHX2
antibody) at dilution of 1:800