

ZMPSTE24 antibody

Product Information

Catalog No.:	FNab09649
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Proteolytically removes the C-terminal three residues of farnesylated proteins. Acts on lamin A/C.

Immunogen information

Immunogen:	zinc metallopeptidase(STE24 homolog, <i>S. cerevisiae</i>)
Synonyms:	CAAX prenyl protease 1 homolog Farnesylated proteins-converting enzyme 1 (FACE-1) Prenyl protein-specific endoprotease 1 Zinc metalloproteinase Ste24 homolog ZMPSTE24 FACE1 STE24
Observed MW:	50 kDa
Uniprot ID :	O75844

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IF: 1:10-1:100
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09649(ZMPSTE24 antibody) at dilution of 1:400