

ZMYND10 antibody

Product Information

Catalog No.:	FNab09650
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for motile ciliary function. Probably involved in axonemal assembly of inner and outer dynein arms(IDA and ODA, respectively) for proper axoneme building for cilia motility. May act by indirectly regulating transcription of dynein proteins.

Immunogen information

Immunogen:	zinc finger, MYND-type containing 10
Synonyms:	Zinc finger MYND domain-containing protein 10 Protein BLu ZMYND10 BLU
Observed MW:	50 kDa
Uniprot ID :	O75800

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000; IP: 1:200-1:1000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab09650(ZMYND10 antibody) at dilution of 1:400