

ZNF124 antibody

Product Information

Catalog No.:	FNab09655
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information


Immunogen:	zinc finger protein 124
Synonyms:	Zinc finger protein 124 Zinc finger protein HZF-16 ZNF124
Observed MW:	40-45 kDa
Uniprot ID :	Q15973

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:10000; IP: 1:1000-1:10000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab09655(ZNF124 Antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab09655(ZNF124 Antibody) at dilution of 1:2000