

ZNF140 antibody

Product Information

Catalog No.:	FNab09657
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation as a repressor.

Immunogen information


Immunogen:	zinc finger protein 140
Synonyms:	Zinc finger protein 140 ZNF140
Observed MW:	53 kDa
Uniprot ID :	P52738

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100
Image:	



Immunofluorescent analysis of HepG2 cells using FNab09657 (ZNF140 antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab09657(ZNF140 Antibody) at dilution of 1:1000