

ZNF174 antibody

Product Information

Catalog No.:	FNab09661
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcriptional repressor.

Immunogen information


Immunogen:	zinc finger protein 174
Synonyms:	Zinc finger protein 174 AW-1 Zinc finger and SCAN domain-containing protein 8 ZNF174 ZSCAN8
Observed MW:	46 kDa
Uniprot ID :	Q15697

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab09661(ZNF174 antibody) at dilution of 1:100



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab09661(ZNF174 antibody) at dilution of 1:500