

## ZNF195 antibody

### Product Information

Catalog No.:	FNab09665
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in transcriptional regulation.

### Immunogen information

Immunogen:	zinc finger protein 195
Synonyms:	Zinc finger protein 195 ZNF195 ZNFP104
Observed MW:	64 kDa
Uniprot ID :	O14628

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab09665(ZNF195 Antibody) at dilution of 1:1000