

## ZNF253 antibody

### Product Information

Catalog No.:	FNab09674
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May function as a transcription factor. Seem to have a transcriptional repression activity.

### Immunogen information

Immunogen:	zinc finger protein 253
Synonyms:	Zinc finger protein 253 Bone marrow zinc finger 1 (BMZF-1) Zinc finger protein 411 ZNF253 BMZF1 ZNF411
Observed MW:	53 kDa
Uniprot ID :	O75346

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09674(ZNF253 Antibody) at dilution of 1:300