

## ZPR1 antibody

### Product Information

Catalog No.:	FNab09676
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


ZNF259, also named as Zpr1, is an evolutionarily conserved protein characterized by two C4 zinc fingers and two Zpr1-specific conserved homology domains. Zpr1 was first identified in mammalian cells as a cytoplasmic protein that was capable of binding the intracellular domain of the unactivated EGF receptor(EGFR). The MW of this protein is 51 kDa, and this antibody specially recognises the 51 kDa protein.

### Immunogen information


Immunogen:	zinc finger protein 259
Synonyms:	Zinc finger protein ZPR1 Zinc finger protein 259 ZPR1 ZNF259
Observed MW:	51 kDa
Uniprot ID :	O75312

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor slide using FNab09676(ZNF259 Antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09676(ZNF259 antibody) at dilution of 1:500