

ZNF266 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab09677 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


May be involved in transcriptional regulation.

Immunogen information


| | |
|--------------|--|
| Immunogen: | zinc finger protein 266 |
| Synonyms: | Zinc finger protein 266 Zinc finger protein HZF1 ZNF266 KIAA2007 |
| Observed MW: | 55 kDa and 70 kDa |
| Uniprot ID : | Q14584 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500-1:5000; IHC: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab09677(ZNF266 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09677(ZNF266 antibody) at dilution of 1:500