

ZNF311 antibody

Product Information

Catalog No.:	FNab09683
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information


Immunogen:	zinc finger protein 311
Synonyms:	Zinc finger protein 311 Zinc finger protein zfp-31 ZNF311 ZFP31
Observed MW:	70 kDa and 60 kDa
Uniprot ID :	Q5JNZ3

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis slide using FNab09683(ZNF311 Antibody) at dilution of 1:50



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab09683(ZNF311 antibody) at dilution of 1:4000