

## ZNF366 antibody

### Product Information

Catalog No.:	FNab09691
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has transcriptional repression activity. Acts as corepressor of ESR1; the function seems to involve CTBP1 and histone deacetylases.

### Immunogen information

Immunogen:	zinc finger protein 366
Synonyms:	Zinc finger protein 366 Dendritic cell-specific transcript protein (DC-SCRIPT) ZNF366
Observed MW:	85 kDa
Uniprot ID :	Q8N895

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:1000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab09691(ZNF366 Antibody) at dilution of 1:300