

ZNF391 antibody

Product Information

Catalog No.:	FNab09693
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 391
Synonyms:	Zinc finger protein 391 ZNF391
Observed MW:	35 kDa
Uniprot ID :	Q9UJN7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:200-1:1000
Image:	



IP Result of anti-ZNF391 (IP:FNab09693, 3ug;
Detection:FNab09693 1:300) with mouse testis
tissue lysate 4000ug.