

## ZNF397 antibody

### Product Information

Catalog No.:	FNab09696
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Isoform 3 acts as a DNA-dependent transcriptional repressor.

### Immunogen information

Immunogen:	zinc finger protein 397
Synonyms:	Zinc finger protein 397 Zinc finger and SCAN domain-containing protein 15 Zinc finger protein 47 ZNF397 ZNF47 ZSCAN15
Observed MW:	61 kDa
Uniprot ID :	Q8NF99

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09696(ZNF397 Antibody) at dilution of 1:300