

ZNF443 antibody

Product Information

Catalog No.:	FNab09699
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 443
Synonyms:	Zinc finger protein 443 Krueppel-type zinc finger protein ZK1 ZNF443
Observed MW:	72 kDa
Uniprot ID :	Q9Y2A4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab09699(ZNF443 antibody) at dilution of 1:500