

## ZNF470 antibody

### Product Information

Catalog No.:	FNab09702
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in transcriptional regulation.

### Immunogen information

Immunogen:	zinc finger protein 470
Synonyms:	Zinc finger protein 470 Chondrogenesis zinc finger protein 1 (CZF-1) ZNF470 CZF1
Observed MW:	69 kDa
Uniprot ID :	Q6ECI4

### Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



rat testis tissue were subjected to SDS PAGE followed by western blot with FNab09702(ZNF470 Antibody) at dilution of 1:1000