

## ZNF483 antibody

### Product Information

Catalog No.:	FNab09704
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in transcriptional regulation.

### Immunogen information

Immunogen:	zinc finger protein 483
Synonyms:	Zinc finger protein 483 Zinc finger protein with KRAB and SCAN domains 16 ZNF483 KIAA1962 ZKSCAN16
Observed MW:	100 kDa
Uniprot ID :	Q8TF39

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab09704(ZNF483 Antibody) at dilution of 1:300