

ZNF496 antibody

Product Information

Catalog No.:	FNab09707
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


DNA-binding transcription factor that can both act as an activator and a repressor.

Immunogen information

Immunogen:	zinc finger protein 496
Synonyms:	Zinc finger protein 496 Zinc finger protein with KRAB and SCAN domains 17 ZNF496 ZKSCAN17
Observed MW:	67 kDa
Uniprot ID :	Q96IT1

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab09707(ZNF496 antibody) at dilution of 1:1000