

ZBTB49 antibody

Product Information

Catalog No.:	FNab09710
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	zinc finger protein 509
Synonyms:	Zinc finger and BTB domain-containing protein 49 Zinc finger protein 509 ZBTB49 ZNF509
Observed MW:	85 kDa
Uniprot ID :	Q6ZSB9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab09710(ZNF509 Antibody) at dilution of 1:600