

## ZNF549 antibody

### Product Information

Catalog No.:	FNab09713
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in transcriptional regulation.

### Immunogen information

Immunogen:	zinc finger protein 549
Synonyms:	Zinc finger protein 549 ZNF549
Observed MW:	74 kDa
Uniprot ID :	Q6P9A3

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:20-1:200
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab09713(ZNF549 Antibody) at dilution of 1:600