

## ZNF557 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab09714  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


May be involved in transcriptional regulation.

### Immunogen information

|              |                                |
|--------------|--------------------------------|
| Immunogen:   | zinc finger protein 557        |
| Synonyms:    | Zinc finger protein 557 ZNF557 |
| Observed MW: | 49 kDa                         |
| Uniprot ID : | Q8N988                         |

### Application

|                       |                                  |
|-----------------------|----------------------------------|
| Reactivity:           | Human                            |
| Tested Application:   | ELISA, WB, IF                    |
| Recommended dilution: | WB: 1:200-1:2000; IF: 1:10-1:100 |
| Image:                |                                  |



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09714(ZNF557 Antibody) at dilution of 1:500