

ZNF593 antibody

Product Information

Catalog No.:	FNab09716
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Negatively modulates the DNA binding activity of Oct-2 and therefore its transcriptional regulatory activity. Could act either by binding to DNA octamer or by interacting with Oct-2. May also be a modulator of other octamer-binding proteins.

Immunogen information

Immunogen:	zinc finger protein 593
Synonyms:	Zinc finger protein 593 Zinc finger protein T86 ZNF593 ZT86
Observed MW:	15 kDa
Uniprot ID :	O00488

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab09716(ZNF593 antibody) at dilution of 1:200