

## ZNF622 antibody

### Product Information

Catalog No.:	FNab09721
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May behave as an activator of the bound transcription factor, MYBL2, and be involved in embryonic development.

### Immunogen information


Immunogen:	zinc finger protein 622
Synonyms:	Cytoplasmic 60S subunit biogenesis factor ZNF622 Zinc finger protein 622 Zinc finger-like protein 9 ZNF622 ZPR9
Observed MW:	54 kDa
Uniprot ID :	Q969S3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab09721(ZNF622 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09721(ZNF622 antibody) at dilution of 1:500