

## ZNF71 antibody

### Product Information

Catalog No.:	FNab09732
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in transcriptional regulation.

### Immunogen information

Immunogen:	zinc finger protein 71
Synonyms:	Endothelial zinc finger protein induced by tumor necrosis factor alpha Zinc finger protein 71 ZNF71 EZFIT
Observed MW:	55 kDa
Uniprot ID :	Q9NQZ8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09732(ZNF71 Antibody) at dilution of 1:1000