

ZNF740 antibody

Product Information

Catalog No.:	FNab09735
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information


Immunogen:	zinc finger protein 740
Synonyms:	Zinc finger protein 740 OriLyt TD-element-binding protein 7 ZNF740 TB7
Observed MW:	22 kDa
Uniprot ID :	Q8NDX6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-ZNF740 (IP:FNab09735, 3ug;
Detection:FNab09735 1:600) with mouse brain
tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with
FNab09735(ZNF740 Antibody) at dilution of
1:600