

ZNHIT1 antibody

Product Information

Catalog No.:	FNab09744
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Seems to play a role in p53-mediated apoptosis induction. Binds to NR1D2 and relieves it of its inhibitory effect on the transcription of APOC3 without affecting its DNA-binding activity.

Immunogen information

Immunogen:	zinc finger, HIT type 1
Synonyms:	Zinc finger HIT domain-containing protein 1 Cyclin-G1-binding protein 1 Zinc finger protein subfamily 4A member 1 p18 Hamlet ZNHIT1 CGBP1 ZNFN4A1
Observed MW:	18 kDa
Uniprot ID :	O43257

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09744(ZNHIT1 antibody) at dilution of 1:1000