

TJP2 antibody

Product Information

Catalog No.:	FNab09752
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	tight junction protein 2(zona occludens 2)
Synonyms:	Tight junction protein ZO-2 Tight junction protein 2 Zona occludens protein 2 Zonula occludens protein 2 TJP2 X104 ZO2
Observed MW:	150 kDa
Uniprot ID :	Q9UDY2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09752(ZO-2 antibody) at dilution of 1:1000