

ZSCAN20 antibody

Product Information

Catalog No.:	FNab09757
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information


Immunogen:	zinc finger and SCAN domain containing 20
Synonyms:	Zinc finger and SCAN domain-containing protein 20 Zinc finger protein 31 Zinc finger protein 360 Zinc finger protein KOX29 ZSCAN20 KOX29 ZNF31 ZNF360
Observed MW:	110 kDa
Uniprot ID :	P17040

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab09757(ZSCAN20 Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09757(ZSCAN20 antibody) at dilution of 1:1500