

## ZSCAN21 antibody

### Product Information

Catalog No.:	FNab09758
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Strong transcriptional activator.

### Immunogen information

Immunogen:	zinc finger and SCAN domain containing 21
Synonyms:	Zinc finger and SCAN domain-containing protein 21 Renal carcinoma antigen NY-REN-21 Zinc finger protein 38 homolog (Zfp-38) ZSCAN21 ZFP38 ZNF38
Observed MW:	54 kDa
Uniprot ID :	Q9Y5A6

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09758(ZSCAN21 antibody) at dilution of 1:1000