

ZWINT antibody

Product Information

Catalog No.:	FNab09762
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Part of the MIS12 complex, which is required for kinetochore formation and spindle checkpoint activity. Required to target ZW10 to the kinetochore at prometaphase.

Immunogen information

Immunogen:	ZW10 interactor
Synonyms:	ZW10 interactor ZW10-interacting protein 1 (Zwint-1) ZWINT
Observed MW:	31 kDa
Uniprot ID :	O95229

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	



human skin tissue were subjected to SDS PAGE followed by western blot with FNab09762(ZWINT antibody) at dilution of 1:300