

ZXDC antibody

Product Information

Catalog No.:	FNab09763
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Cooperates with CIITA to promote transcription of MHC class I and MHC class II genes.

Immunogen information

Immunogen:	ZXD family zinc finger C
Synonyms:	Zinc finger protein ZXDC ZXD-like zinc finger protein ZXDC ZXML
Observed MW:	75 kDa
Uniprot ID :	Q2QGD7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab09763(ZXDC antibody) at dilution of 1:300