

---

## Goat Anti-Rabbit IgG(H+L)

### Product Information

Catalog No.:	FNab09778
Size:	1mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Rabbit IgG
Observed MW:	25 kda, 50 kDa

### Application

Reactivity:	rabbit
Tested Application:	ELISA, WB