

STAG3 antibody

Product Information

Catalog No.:	FNab09779
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	stromal antigen 3
Synonyms:	Cohesin subunit SA-3 SCC3 homolog 3 Stromal antigen 3 Stromalin-3 STAG3
Observed MW:	155 kDa
Uniprot ID :	Q9UJ98

Application

Reactivity:	Human
Tested Application:	ELISA, IP, IHC, IF
Recommended dilution:	IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab09779(STAG3 Antibody) at dilution of 1:200



IP Result of anti-STAG3 (IP: FNab09779, 4ug; Detection: FNab09779 1:300) with human testis tissue lysate 4000ug.