

PCDH10 antibody

Product Information

Catalog No.:	FNab09785
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential calcium-dependent cell-adhesion protein.

Immunogen information


Immunogen:	protocadherin 10
Synonyms:	Protocadherin-10 PCDH10 KIAA1400
Observed MW:	113 kDa
Uniprot ID :	Q9P2E7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC,IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab09785(PCDH10 Antibody) at dilution of 1:200



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09785(PCDH10 Antibody) at dilution of 1:300