

FTO antibody

Product Information

Catalog No.:	FNab09787
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None


Immunogen information

Immunogen:	fat mass and obesity associated protein
Synonyms:	Alpha-ketoglutarate-dependent dioxygenase FTO Fat mass and obesity-associated protein U6 small nuclear RNA (2'-O-methyladenosine-N(6)-)-demethylase FTO U6 small nuclear RNA N(6)-methyladenosine-demethylase FTO mRNA (2'-O-methyladenosine-N(6)-)-demethylase FTO (m6A(m)-demethylase FTO) mRNA N(6)-methyladenosine demethylase FTO tRNA N1-methyl adenine demethylase FTO FTO KIAA1752
Observed MW:	58 kDa
Uniprot ID :	Q9C0B1


Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:500

Image:



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab09787(FTO Antibody) at dilution of 1:400



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09787(FTO Antibody) at dilution of 1:300