

## ITIH4 antibody

### Product Information


Catalog No.:	FNab09796
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	ITI heavy chain H4
Synonyms:	Inter-alpha-trypsin inhibitor heavy chain H4 (ITI heavy chain H4, ITI-HC4, Inter-alpha-inhibitor heavy chain 4) Inter-alpha-trypsin inhibitor family heavy chain-related protein (IHRP) Plasma kallikrein sensitive glycoprotein 120 (Gp120, PK-120) 70 kDa inter-alpha-trypsin inhibitor heavy chain H4 35 kDa inter-alpha-trypsin inhibitor heavy chain H4 ITIH4 IHRP ITIHL1 PK120
Observed MW:	120 kDa
Uniprot ID :	Q14624

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using FNab09796(ITIH4 Antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab09796(ITIH4 Antibody) at dilution of 1:600