

FTH1 antibody

Product Information

Catalog No.:	FNab09810
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4A9
IsoType:	IgG1
Storage:	20mMPBS, pH7.4 with 0.1% proclin300, 2-8°C for 2 weeks, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	Ferritin
Synonyms:	Ferritin heavy chain (Ferritin H subunit) Cell proliferation-inducing gene 15 protein Ferritin heavy chain, N-terminally processed FTH1 FTH FTHL6
Observed MW:	170 kDa
Uniprot ID :	P02794

Application

Reactivity:	Human
Tested Application:	ELISA, CLIA, GICA(capture antibody)