

---

## anti- TLR2 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab09836   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Protein A+G purification  |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Mouse   |
| Clonality:    | monoclonal  |
| Clone ID:     | 3H6   |
| IsoType:      | IgG2a   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | toll-like receptor 2                    |
| Synonyms:    | CD282, TIL4, TLR2, toll like receptor 2 |
| Observed MW: | 90 kDa                                  |
| Uniprot ID : | O60603                                  |

### Application

|                          |              |
|--------------------------|--------------|
| Reactivity:              | Human        |
| Tested Application:      | ELISA, WB    |
| Recommended dilution:WB: | 1:500-1:2000 |
| Image:                   |              |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



THP-1 cells were subjected to SDS PAGE followed by western blot with FNab09836 (TLR2 antibody) at dilution of 1:2000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)