
Goat Anti- Mouse IgG(H+L)

Product Information

Catalog No.:	FNab09848
Size:	1mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	mouse IgG
Observed MW:	155 kDa
Uniprot ID :	None

Application

Reactivity:	Mouse
Tested Application:	ELISA, WB, IHC