

RAPGEF3 antibody

Product Information

Catalog No.:	FNab09888
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Rap guanine nucleotide exchange factor(GEF) 3
Synonyms:	Rap guanine nucleotide exchange factor 3 Exchange factor directly activated by cAMP 1 Exchange protein directly activated by cAMP 1 (EPAC 1) Rap1 guanine-nucleotide-exchange factor directly activated by cAMP cAMP-regulated guanine nucleotide exchange factor I (cAMP-GEFI) RAPGEF3 CGEF1 EPAC EPAC1
Observed MW:	104 kDa
Uniprot ID :	O95398

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500 - 1:2000; IF: 1:50 - 1:200
Image:	



Immunofluorescence analysis of U2OS cells using FNab09888(EPAC1 antibody) at dilution of 1:100. Blue: DAPI for nuclear staining.



mouse heart were subjected to SDS PAGE followed by western blot with FNab09888(EPAC1 antibody) at dilution of 1:1000